

of the Atomic Absorption Methods section of this manual.

6.0 Notes

- 6.1 High levels of silicon may interfere
- 6.2 The air-acetylene flame absorbs about 25% of the energy at the 213.9 nm line.
- 6.3 The sensitivity may be increased by the use of low-temperature flames.
- 6.4 Some sample container cap liners can be a source of zinc contamination. To circumvent or avoid this problem, the use of polypropylene caps is recommended.
- 6.5 The dithizone colorimetric method may also be used (Standard Methods, 14th Edition, p 265).
- 6.6 For concentrations of zinc below 0.01 mg/L, either the Special Extraction Procedure given in part 9.2 of the Atomic Absorption Methods section or the furnace procedure, Method 289.2, is recommended.
- 6.7 Data to entered into Storet must be reported as $\mu\text{g/L}$.

7.0 Precision and Accuracy

- 7.1 An interlaboratory study on trace metal analyses by atomic absorption was conducted by the Quality Assurance and Laboratory Evaluation Branch of EMSL. Six synthetic concentrates containing varying levels of aluminum, cadmium, chromium, copper, iron, manganese, lead and zinc were added to natural water samples. The statistical results for zinc were as follows:

Number of Labs	True Values $\mu\text{g/Liter}$	Standard Mean Value $\mu\text{g/Liter}$	Deviation $\mu\text{g/Liter}$	Accuracy as % Bias
86	281	284	97	1.2
89	310	308	114	-0.7
82	56	62	28	11.3
81	70	75	28	6.6
62	7	22	26	206
61	11	17	18	56.6