



- 5.3 Dual-Purpose Sulfite Indicator Powder: a proprietary formulation containing sulfamic acid to destroy nitrite.
- 5.4 Standard potassium iodide-iodate titrant, 0.0125N: Dissolve 445.8 mg anhydrous potassium iodate,  $\text{KIO}_3$  (primary standard grade dried for several hours at  $120^\circ\text{C}$ ), 4.25 g KI and 310 mg  $\text{NaHCO}_3$  in distilled water and dilute to 1 liter. This titrant is equivalent to  $500 \mu\text{g SO}_3/1.00 \text{ mL}$ .
- 5.5 Sulfamic Acid: Crystalline
- 5.6 EDTA Reagent Dissolve 2.5 g EDTA in 100 mL distilled water.

## 6.0 Procedure

- 6.1 Sampling Contact with air must be minimized. If the sample temperature is greater than  $50^\circ\text{C}$ , it must be cooled in a special apparatus described elsewhere (see Bibliography). Immediately add 1 mL of EDTA Solution (5.6) per 100 mL of sample.
- 6.2 Starch Indicator
  - 6.2.1 Place 1 mL  $\text{H}_2\text{SO}_4$  (5.1) in titration vessel.
  - 6.2.2 Add 0.1 g sulfamic acid crystals (5.5).
  - 6.2.3 Add 50 mL sample.
  - 6.2.4 Add approximately 0.1 g starch indicator (5.2).
  - 6.2.5 Titrate with potassium iodide-iodate titrant (5.4) until a faint permanent blue color develops. Keep the pipet tip below the surface of the sample. View the color change against a white background. Record the mL titrant.
  - 6.2.6 Run a reagent blank using distilled water instead of sample (6.2.3).
- 6.3 Dual Purpose Sulfite Indicator Powder.
  - 6.3.1 Place 50 mL sample in a titration vessel.
  - 6.3.2 Add 3-4 drops phenolphthalein indicator.
  - 6.3.3 Add sufficient scoops (1 g) of indicator (5.3) to discharge the red color.
  - 6.3.4 Titrate with potassium iodide-iodate titrant (5.4) until a faint permanent blue color develops. View the color change against a white background. Record the mL titrant.
  - 6.3.5 Run a reagent blank using distilled water instead of sample (6.3.1).

## 7.0 Calculations

- 7.1 Use the formula:

$$\text{mg/L SO}_3 = \frac{A \times N \times 40,000}{\text{mL sample}}$$

where:

A = mL titrant (6.2.5 or 6.3.4) B = mL  
 B = mL titrant for the blank (6.2.6 or 6.3.5) and  
 N = normality of KI-KIO<sub>3</sub> titrant (5.4)

- 7.2 To calculate as  $\text{Na}_2\text{SO}_3$   
 $\text{mg/L Na}_2\text{SO}_3 = \text{mg/L SQ} \times 1.57$

## 8.0 Precision and Accuracy

8.1 Precision and accuracy data are not available at this time.

### **Bibliography**

1. Annual Book of ASTM Standards, Part 31, "Water", Standard D1339-72, Method C, p 440 (1976).
2. Standard Methods for the Examination of Water and Wastewater, 14th Edition, p 508, Method 429, (1975).